



RenseWAI GE

RenseWAI GE is an EDI (Electro De-Ionization) based Pure Water System, produces Type 2 or EDI water directly from potable tap water. Resistivity of the product water is above 5 M Ω •cm at 25°C, which meets or exceeds Type II water quality, as defined by ASTM reagent grade water standards.

Features

Design & Functions

- The system can be linked to multiple dispensers via CAN cable or wireless.
- Automatic temperature compensation allows RO production rate stable over a wide range of temperature.
- Exceptionally consistent and predictable high purity Type II water, from the best in class IonPure EDI (Electro De-Ionization) module.
- Fully automated control system has cleaning, rising, flush and many other functions.
- System automatically rinses a new RO membrane.
- A bypath automatically sends the RO water to the drain if it does not meet a preset quality standard. The RO drain water is reused / recycled to increase the water yield thus it is more environmentally friendly.
- The main system can be either wall-mounted or set on the shelf or under the sink, to save precious lab space.
- The Chlorine cleaning and pH cleaning functions keep the maintenance easier and ensure optimal water quality.





Effective TOC level control

- Built-in 254 nm lamp kills bacteria in pure water.
- A tank Sanitization Module can be equipped to effectively keep microbial growth at a lower level by UV light. (optional)

Easy Operation and Maintenance

- Performance history, alarms and alerts information, maintenance data of consumables as well as key parts are logged and traceable from a simple RFID scan at installation.
- RFID tags ensure perfect placement of consumables and trace their performance.
- Water quality, operation parameters, the status of the system, dispensers, components, and peripheral devices are stored and displayed on the large color touch screen monitor.
- Signature verification for maintenance and service.

Cartridges and Parts

1. EDI module

- Internal Ion-Pure EDI module removes the remaining ions, further reduces the content of organic compounds.
- Ion exchange resins are continuously regenerated by an electric field. No hazardous chemical regeneration or costly resin replacement is needed.
- EDI module has a high degree of automation and minimal electricity consumption by using the equivalent of the energy.

2. Cartridge

- Internal P Pack cartridge removes organics, particles and scaling ions to prevent them from fouling the RO membranes and the EDI module.
- AC Pack contains special activated carbon to remove organic impurities and free chlorine.

3. Control console

- An 8-inch console controls system and peripheral devices. All operations can be done on the console by touching the screen.
- The control console can be placed on the bench or inside a drawer for further flexibility and space saving.
- The console screen and dispenser screens are water-proof. You can operate the console and dispenser with latex gloves on.
- The system has two level password protection on display menu for safety & security purposes.

4. Dispenser

- The volumetric function can automatically deliver the water volume you need once it is set up.
- With the adjusting button on a control console or a dispenser, you can dispense water at the flow rate you need up to 2 liters per minute.
- A dispenser handle can be set on the dispenser stand, on the main system, or even hanging onto other places to free up maximum bench space.





5. Final Filter

• Various final filters (optional) are applicable to ensure ultrapure water without particles, bacteria or pyrogen.

6. Level sensor

- The tank integrated continuous level sensor measures water level within the tank and manages the system to start or stop producing water automatically based on requirements.
- Water level and quality inside the tank can be viewed from the control console.

Main Applications

EDI Pure Water

- Preparation of chemical bio-reagents
- Preparation of culture media
- Preparation of solution for chemical analysis such as HPLC and ICP
- Feed water to ultrapure water systems
- Feed water to medical device and equipment (clinical analyzers, Aging testers etc.)
- For serum and blood fractionation
- For ophthalmic

Main Specification

Feed Water Requirements

- ➤ Feed Water Potable Tap Water
- \triangleright Feed water conductivity < 2000 µS/cm or TDS < 1000 ppm
- ➤ Feed water pressure (2-6) bar
- ➤ Operating temperature (5 35)°C
- > Feed water TOC < 1000 ppb (preferable)

Flow rate

➤ Type 2 or EDI based production rate – 5, 10 and 15 L/hr.

Product Water Quality

- > Type 2 or EDI based water resistivity (@25°C) > 5 MΩ·cm (typically 10 16 MΩ·cm)
- > Type 2 or EDI based water TOC < 30 ppb


Power

- Input Voltage 110 240 VAC
- Operating Voltage 24 VDC
- ➤ Wattage < 200 W





Flow Chart



Ordering Information

SI.No.	Description	Catalogue No.
1	RenseWAI GE 5 System SET	RW0E00500K
2	RenseWAI GE 10 System SET	RW0E01000K
3	RenseWAI GE 15 System SET	RW0E01500K





Standard contents of the Set

RenseWAI GE System Set consists of :-

- 1. Main system
- 2. Control Console
- 3. RO membrane
- 4. UV lamp 254 nm (pre-installed)
- 5. AC Pack
- 6. P Pack
- 7. Reservoir 30 litre capacity, with tank level sensor
- 8. Tank vent filter (with Soda lime in granular form, which removes carbon dioxide CO2 and Granular activated carbon or GAC, which removes organic matter from the air and 0.22 hydrophobic membrane, which removes air-borne bacteria)
- 9. Three stage prefiltration kit (PF Kit)
- (1 micron + 10 micron + Carbon Cartridge of 3 micron pore size, included inside PF Kit.Quantity -1 no. each, all 10'' long)
- 10. External feed booster pump, with high & low pressure auto cutoff switch included.

(Included / needed only in-case, if potable feed water pressure is less than 2 bar)

Photographs (Representative only)












(Disclaimer – Please note photographs are for representation purpose only & can vary from the actual system)



RenseWAI is the registered trademark of M/s. Pramuk Healthcare



Rensewa





Rensewa











